

Vacutec™ EV

- Low level noise at less than 58 dBA.
- Integrated base reduces vibration.
- Designed for single patient, everyday use.
- Easy-to-read gauge and easy-to-turn flow regulator knob.
- Suction cup feet for stable operation.
- MET certified and complies with UL and CSA listings.
- Complete with disposable 800cc collection bottle.
- Automatic shut-off.
- Designed to meet the requirements of HPSC code: E0600
- Two years limited warranty.

JB0112-014 1 ea

ITEM	Width	Depth	Height	Weight	Power Requirements	Vacuum Range
JB0112-014	8 3/4"	14 1/2"	10 1/4"	12.1 lb	115 VAC, 60 Hz	0 – 560 mmHg

